

New GeneiPureID™ DNA Isolation Kits


Description: The ability to prepare and isolate pure nucleic acids from a variety of sources is an important step in many molecular biology protocols. The GeneiPureID™ DNA Isolation Kit provides a simple, fast, reliable, reproducible and convenient extraction of high amount of pure genomic DNA from various forensic samples (saliva, buccal swab, hair, nail and dried blood). The kit uses specialized lysis solution (specific for the different tissues) and spin columns for purifying the DNA thus obtained.

Protocol for GeneiPureID™ DNA IsolationSample processing
↓


Lysis of samples



Bind



Wash



Elute