

## Instruments for PCR based Detection Kits

Bangalore Genei can supply laboratory equipment to start a new PCR laboratory, for WSSV, STWSSV & MBV detections. Please refer to the ordering information.

### Ordering Information

Catalogue Number	Product	Pack Size
<b>INSTRUMENTS FOR PCR BASED DETECTION KITS</b>		
107214	Midi PCR Workstation (HxWxD: 45x60x30) cm (Designed for UV Irradiation of PCR Chemicals & Instruments)	1 Set
106888	Miniphor UVT System Complete Set with 8, 12, 16 well combs	1 Set
106900	Cooling unit for Miniphor UVT System	1 Set
106856	<b>GeNei™</b> Digital PS 100: EPS 300 Model ; 0-300 V; 0-100 mA ; Variable with timer Constant Current ; Constant Voltage ; Digital display	1 Unit
107161	UV Transilluminator with filter (312 nm) 20 x 20 cm Imported	1 No.
107145	UV Face Shield Imported	1 No.
107108	Handheld UV Lamp 6 Watt 50 x150 mm (365, 254 nm) Imported	1 No.
107100	<b>GeNei™</b> Cooler (24 x 1.5 ml tube; An alternate to ice buckets)	1 Unit
107102	<b>GeNei™</b> Cooler (24 x 0.5 ml tube; An alternate to ice buckets)	1 Unit
107944	Microcentrifuge (fixed rotor speed 10000 rpm) with digital timer	1 Unit
107170	Dry bath with standard heating block	1 Set
106915	Micropipette 2-20 $\mu$ l increment 0.1 $\mu$ l	1 No.
106975	Universal fit aerosol filter barrier tips, volume 20 $\mu$ l	1000/pack
106930	Micropipette 20-200 $\mu$ l increment 1.0 $\mu$ l	1 No.
106937	0-200 $\mu$ l Clear tips	1000/pack
106931	0.2 ml Thin Wall PCR Tubes - Domed	1000 tubes
106932	0.5 ml Thin Wall PCR Tubes, Flat Cap	1000 tubes
106933	Storage rack for PCR-0.2 ml (96 wells)	1 No.