

## New GeneiPure™ Bacterial DNA Purification Kit

**Highlights:**

- Improved yield, upto 8-10 µg of DNA from 0.5-1 ml of *Escherichia coli* culture
- Unique lysis buffer optimized for Bacterial DNA Purification

The kit provides sufficient reagents for 20/50 preparations.

**Storage:** Store cell lysis solutions at 4°C, enzymes at -20°C and rest of the components at RT.


Material	Amount	Yield (µg)
<i>E.coli</i>	1.5 ml	5-7
<i>Bacillus subtilis</i>	4.5 ml	10-13

**Materials provided:**

1. Cell Lysis Solutions
2. Lysozyme
3. Proteinase K
4. RNase A
5. Wash Buffers
6. Elution Buffer
7. GeneiPure™ Columns
8. 2 ml Collection Tubes
9. Instruction Manual

**Ordering Information:**

Product	Size	Cat #
GeneiPure™ Bacterial DNA Purification Kit (20 Preparations)	1 Pack	117290
GeneiPure™ Bacterial DNA Purification Kit (50 Preparations)	1 Pack	117291



Lane M : Lambda/Hind III Marker  
 Lane 1 : *Staphylococcus aureus* 3A  
 Lane 2 : *Bacillus amyloliquefaciens*  
 Lane 3 : *Bacillus globigi*  
 Lane 4 : *Klebsiella pneumonia*  
 Lane 5 : *Neisseria sicca*  
 Lane 6 : *Escherichia coli* C