

## New GeneiPure™ Yeast DNA Purification Kit

### Highlights:

- Improved yield, upto 5-8 µg of DNA from 4-5 ml of *Saccharomyces cerevisiae* culture
- Unique lysis buffer and lyticase for effective and complete lysis

The kit provides sufficient reagents for 20/50 preparations.

**Storage:** Store cell lysis solutions at 4 °C, enzymes at -20 °C and rest of the components at RT.

Material	Amount	Yield (µg)
<i>Saccharomyces cerevisiae</i>	4-5 ml	5-8

### Ordering Information:

Product	Size	Cat #
GeneiPure™ Yeast DNA Purification Kit (20 Preparations)	1 Pack	117295
GeneiPure™ Yeast DNA Purification Kit (50 Preparations)	1 Pack	117294

### Materials provided:

1. Cell Lysis Solutions
2. Lyticase
3. Proteinase K
4. RNase A
5. Wash Buffers
6. Elution Buffer
7. GeneiPure™ Columns
8. 2 ml Collection Tubes
9. Instruction Manual