

New GeneiPure™ Plant Genomic DNA Purification Kit

Highlights:

- Upto 100 mg of plant sample can be processed with a typical yield of 10 µg of DNA
- Two alternate lysis buffers for processing of different plant species
- Improved lysis buffers for high yield and purity
- Spin filters included for clarification of lysate

The kit provides sufficient reagents for 20/50 preparations.

Storage: Store lysis solutions at 4 °C, enzymes at -20 °C and rest of the components at RT.

Materials provided:

1. Cell Lysis Solutions
2. Elution Buffer
3. Proteinase K
4. RNase A
5. GeneiPure™ Columns (20 Nos.)
6. Wash Buffer
7. Tissue Grinder
8. Instruction Manual

Ordering Information:

Product	Size	Cat #
GeneiPure™ Plant Genomic DNA Purification Kit (20 Preparations)	1 Pack	117298
GeneiPure™ Plant Genomic DNA Purification Kit (50 Preparations)	1 Pack	117297