

New GeneiUltrapure™ Plasmid Mini & Midi Purification Kit

Highlights:

- Fast & Economical purification
- Available in two pack sizes Mini & Midi
- Purification of 2-25kb plasmids
- Recovery >90%
- Suitable for both High and Low copy plasmids
- Special folded filters supplied only with Midi Kit
 - Helps removal of Cell debris
 - Prevents shearing of plasmid DNA
 - Replaces time consuming centrifugation step

Materials Provided:

1. Buffers for Alkali lysis (Solution I, II & III)
2. Buffers for Columns (Equilibration, Wash & Elution)
3. RNaseA (Lyophilized)
4. GeneiUltraPure™ Columns
5. Instruction Manual

Storage: All components store at RT.

Ordering Information:

Product	Size	Cat #
GeneiUltraPure™ Plasmid Mini Purification Kit (20 Purifications)	1 Pack	117766
GeneiUltraPure™ Plasmid Midi Purification Kit (10 Purifications)	1 Pack	117767

Specification:

Product	Size	Capacity	Culture volume
GeneiUltrapure™ Plasmid Mini Prep Kit	20 Purifications	20 µg	Low Copy-3 ml
			High Copy-10 ml
GeneiUltrapure™ Plasmid Midi Prep Kit	10 Purifications	100µg	Low copy -25 ml
			High Copy-100 ml