

**New GeneiUltraPure™ Bacterial Genomic DNA Purification Kit**

**Highlights:**

- Suitable for both Gram Positive and Gram Negative Bacteria
- Lysozyme provided for effective lysis

**Ordering Information:**

Product	Size	Cat #
GeneiUltraPure™ Bacterial Genomic DNA Purification Kit (10 Purifications)	1 Pack	117771

**Materials Provided:**

1. Cell Lysis Buffers
2. Buffers for Columns (Equilibration, Wash & Elution)
3. RNaseA (Lyophilized)
4. Lysozyme (Lyophilized)
5. ProteinaseK (Lyophilized)
6. GeneiUltraPure™ Columns
7. Instruction Manual

**Storage:** Cell lysis buffers at 4°C, enzymes at -20°C and rest of the components at RT

**Specification:**

Product	Size	Amount	Yield
GeneiUltraPure™ Bacterial Genomic DNA Purification Kit	10 Purifications	0.5-1 ml (gram-ve)	10-15 µg
		1-2 ml (gram + ve)	10-15 µg