

## New GeneiUltrapure™ Plant Genomic DNA Purification Kit

### Highlights:

- Suitable for wide range of Plant Species

### Ordering Information:

### Materials Provided:

1. Cell Lysis Buffers
2. Buffers for Columns (Equilibration, Wash & Elution)
3. RnaseA (Lyophilized)
4. ProteinaseK (Lyophilized)
5. GeneiUltraPure™ Columns
6. Instruction Manual

Product	Size	Cat #
GeneiUltraPure™ Plant Genomic DNA Purification Kit (10 Purifications)	1 Pack	117772

**Storage:** Cell lysis buffers at 4°C, enzymes at -20°C and rest of the components at RT

### Specification:

Product	Size	Material	Dry tissue (mg)	Yield (in µg)
GeneiUltrapure™ Genomic Prep Kit from Plant	10 Purifications	Tea	100	3-5
		Papaya	100	5-7
		Banana	100	2-3