

New GeneiUltraPure™ Cotton Genomic DNA Purification Kit - Leaves

Highlights:

- High yields from Transgenic or non transgenic Cotton Leaves
- Higher yields obtained for effective lysis using additives

Materials Provided:

1. Cell Lysis Buffers
2. Buffers for Columns (Equilibration, Wash & Elution)
3. RNase A (Lyophilized)
4. Proteinase K (Lyophilized)
5. Additive
6. GeneiUltraPure™ Columns
7. Instruction Manual

Storage: Cell lysis buffers at 4°C, enzymes at -20°C and rest of the components at RT

Ordering Information:

Product	Size	Cat #
GeneiUltraPure™ Cotton Genomic DNA Purification Kit - Leaves (10 Purifications)	1 Pack	117773

Specification:

Product	Size	Amount	Yield
GeneiUltraPure™ Cotton Genomic DNA Purification Kit - Leaves	10 Purifications	100 mg (wet tissue)	3-5 µg