

New GeneiUltraPure™ Mammalian Genomic DNA Purification Kit - Tissue

Highlights:

- Isolation of Genomic DNA from wide range of tissues
- Suitable for both Fresh and frozen Tissues
- Sterilized Tissue Homogenizers provided.

Materials Provided:

1. Cell Lysis Buffers
2. Buffers for Columns (Equilibration, Wash & Elution)
3. RNase A (Lyophilized)
4. Proteinase K (Lyophilized)
5. Tissue Homogenizer
6. GeneiUltraPure™ Columns
7. Instruction Manual

Storage: Cell lysis buffers at 4°C, enzymes at -20°C and rest of the components at RT

Ordering Information:

Product	Size	Cat #
GeneiUltraPure™ Mammalian Genomic DNA Purification Kit - Tissue (10 Purifications)	1 Pack	117769

Specification:

Product	Size	Amount	Yield
GeneiUltraPure™ Mammalian Genomic DNA Purification Kit Tissue	10 Purifications	10-20 mg Tissues	20-30 µg
		2x10 ⁶ cells	20-30 µg