

New GeneiUltrapure™ Blood Genomic DNA Purification Kit

Highlights:

- Isolation of Genomic DNA from whole blood
- Can be used for both fresh and frozen blood
- Complete removal of RBC

Materials Provided:

1. Cell Lysis Buffers
2. Buffers for Columns (Equilibration, Wash & Elution)
3. RNaseA (Lyophilized)
4. ProteinaseK (Lyophilized)
5. GeneiUltraPure™ Columns
6. Instruction Manual

Storage: Cell lysis buffers at 4°C, enzymes at -20°C and rest of the components at RT

Ordering Information:

Product	Size	Cat #
GeneiUltraPure™ Blood Genomic DNA Purification Kit (10 Purifications)	1 Pack	117768

Specification:

Product	Size	Amount	Yield
GeneiUltrapure™ Blood Genomic DNA Purification Kit	10 Purifications	300 µl	2-5 µg