

New Monoclonal IgG Purification Kit

Description: GeNei™ Monoclonal IgG purification kit is based on the affinity of IgG in monoclonal culture supernatant or ascites fluid to bind to Protein A-Sepharose at high salt concentration and elute at low pH. The kit enables easy isolation of highly pure IgG1 fraction from the above sources.

The kit contains all materials and reagents required for purification of 3 batches of IgG.

Materials Provided:

- Protein A- Sepharose Affinity Column
- Equilibration Buffer
- Elution Buffer
- Neutralizing Buffer
- Storage Buffer
- Instruction Manual

Storage: 4°C.

Ordering Information:

Product	Size	Cat #
Monoclonal IgG Purification Kit	1 Pack	117746

GeNei™ Immunoglobulin Purification Kit (Thiophilic Column)

Description: GeNei™ Immunoglobulin purification kit is a unique kit for extracting immunoglobulins from serum of mammalian or avian species. The kit is based on the ability of thiophilic column to bind immunoglobulins from variety of species.

Materials Provided:

- Thiophilic column of (5 ml)
- Equilibration reagent
- Elution reagent
- Instruction Manual.

Storage: 4°C.

Ordering Information:

Product	Size	Cat #
GeNei™ Immunoglobulin Purification Kit (Thiophilic Column)	1 No.	105576