

## Vertical Gel Electrophoresis System(Imported)

**Description:** Vertical Gel Electrophoresis System complete system, system includes electrode module, electrophoresis tank, casting module-2 sets, glass plate set-5 nos, 0.5mm spacer-4 nos, (0.5 mm or 0.75 mm or 1.0 mm) x 10 well comb), One gel gate, Alignment card, glass plate separator, instruction manual.

- Optional Accessories:**
- Electrode Module
  - Casting Module
  - Glass plate Set
  - 1.0mm Spacer
  - 1mm 10 teeth comb
  - Alignment card
  - Glass plate Separator

## Paper Electrophoresis System

**Description:** Basic unit consists of anode and cathode buffer reservoirs capable of accommodating suitable electrode assembly. In case of Vertical model, paper supporting vertical stands with cross rods are provided in between reservoirs, whereas, a removable horizontal supporting stand has been provided in horizontal model. In both the cases transparent acrylic cover are provided for controlling buffer evaporation and for electrical safety.

**Vertical Mini for 3 strips:** Basic unit with lid (1 Set) Platinum Electrodes (2 nos.), Paper Supporting Rods (3 nos), Connecting Cord (1 set).

**Vertical Large for 6 strips:** Vertical large is for 6 Strips 10 Micro slides. The system contains Basic unit with lid (1set), Platinum Electrodes (2 nos), Paper supporting rod (3 nos.) Connecting Cord (1 set).

## Immuno electrophoresis Systems

Electrophoresis Systems for various protein and immuno electrophoresis technique.

### Mini Immuno electrophoresis system (Cat # 106402)

**Description:** Basic unit (1 No.), platinum electrode assembly detachable (2 Nos.), Lid (1 No.), 1-well comb with supporting stand (1 No.), glass plates 50 x 50 mm (6 Nos.), glass plates 75 x 50 mm (6 Nos.), Acrylic template with well former (1 No.), gel puncher 4.0 mm o.d. (1 No.) connecting cord (2 Nos.) and instruction manual.

### Midi Immuno electrophoresis system (Cat # 106413)

**Description:** Basic unit (1 No.), platinum electrode assembly detachable (2 Nos.), glass plates 75 x 50 mm (6 Nos.), glass plates 50 x 50 mm (6 Nos.), CIE template (1 set), gel puncher 4.0 mm (o.d.) (1 No.), Acrylic Template with well former and 'U' spacers (1 set), lid (1 No.), connecting cord (2 Nos.) and instruction manual.

### Immuno/Submarine Multi Purpose Model (Cat # 106419)

**Description:** Basic unit (1 No.), lid (1 No.), platinum electrode (2 Nos.), connecting cord (2 Nos.), 8 well comb, 1.5 mm (2 Nos.), 13 well comb, 1.5 mm (1 No.) comb stand (1 No.), gel platform 130 x 130 mm (1 No.), gel platform 130 x 60 mm (2 Nos.), levelling table (1 No.), glass plates 100 x 80 mm (6 Nos.), gel punchers 3, 4, 5 mm (o.d.) (1 each), clips (4 Nos.), variable radius template (1 No.), 5-well template (1 No.), 8-well template (1 No.), multiple template with holder and spacer and glass plate 200 x 90 mm (1 set), Graber Williams template with holder and spacer and glass plates 100 x 150 mm (1 set), gel cutter (1 No.), gel remover (1 No.), glass plate and template for electrofocusing 130 x 260 mm (1 set).

**Accessories:** Please refer to the ordering information

#### Recommended Power Supply:

##### Digital Power Supply

**GeNei™** Digital PS 500 : 500 V / 500 mA (Cat # 106857)/Programmable Power supply (500 V, 500 mA) (Cat # 116740)

## Ordering Information:

Catalogue Number	Product	Pack Size	Price (Rs.)
<b>VERTICAL GEL ELECTROPHORESIS SYSTEM (IMPORTED)</b>			
116667	Vertical Gel Electrophoresis System (Imported) (with (0.5 mm or 0.75 mm or 1.0 mm) x 10 well comb)	1 Unit.	59500
<b>PAPER ELECTROPHORESIS SYSTEM</b>			
116723	Paper Electrophoresis System (Vertical Mini for 3 strips)	1 Unit.	11250
116721	Paper Electrophoresis System (Vertical Large for 6 strips)	1 Unit.	16100
<b>IMMUNOELECTROPHORESIS SYSTEMS</b>			
106402	Mini Immunoelectrophoresis System	1 Unit	5000
106413	Midi Immunoelectrophoresis System	1 Unit	8350
106419	Immuno / Submarine Multipurpose System	1 Unit	15000
<b>ACCESSORIES for Mini Immunoelectrophoresis System</b>			
107215	1-Well comb with stand	1 Set	700
106410	Template with spacer	1 Set	650
106411	Glass plates (50 mm x 50 mm) (12/pack)	1 Pack	650
106412	Glass plates (75 mm x 50 mm) (12/pack)	1 Pack	700
<b>For Midi and Multipurpose System</b>			
106414	Immuno diffusion template- 75 mm x 50 mm	1 No.	650
106415	CIE template	1 No.	650
106417	Gel punchers (3,4,5 mm o.d)	1 No. each	700
106418	Levelling table with levelling screws and spirit level	1 No.	1800
106437	Electrofocussing plate with 0.5 mm spacer	1 No.	1300
106438	Electrofocussing plate with 1.0 mm spacer	1 No.	1300
106440	Electrofocussing plate with 1.5 mm spacer	1 No.	1300
106527	Gel running tray (130 mm x 130 mm)	1 No.	1000
106556	Gel casting tray (130 mm x 130 mm)	1 No.	1200
106441	Gel running tray (130 mm x 60 mm) (2 / pack)	1 Pack	1000
106447	Gel casting tray (130 mm x 60 mm)	1 No.	1300
106539	Acrylic comb 13 well 1.5 mm (In Built)	1 No.	700
106529	Acrylic comb 8 well 1.5 mm (In Built)	1 No.	650
106442	Acrylic comb preparative 3-well 3.0 mm (In Built)	1 No.	700
106443	Acrylic comb preparative 3-well 6.0 mm (In Built)	1 No.	800
106444	Developing tray for PAGE	1 No.	1500
106445	Developing tray for agarose gel (light-proof)	1 No.	1300
106446	Universal punchers & templates	1 Set	2900
106994	Connecting cord (Red & Black)	1 Set	550
106448	Glass plates (100 mm x 80 mm)	6/Pack	700
106449	Glass plates (100 mm x 115 mm)	6/Pack	800
106460	Glass plates (210 mm x 90 mm)	6/Pack	900
106461	Glass plates (260 mm x 130 mm)	6/Pack	1000