
Pipettes & Pipette Stand

Digital Variable Volume Pipette (Autoclavable)

- These imported pipettes are convenient, easy to use, having locking mechanism that ensures the selected volume, remains securely fixed and does not change during operation.
- The preset volume is confirmed by the convenient digital volume indicator.
- The entire pipette can be autoclaved (20 minutes at 121° C) for decontamination and sterilization.

Digital Variable Volume Pipette - Popular Model

These micropipettes are designed for use in Educational, Research, Diagnostics and Clinical pathology Labs. **The tip cone** of these micropipettes is fully autoclavable which keeps it sterile for each use.

Pipette Stands

- Pipette stand can hold 5 standard size pipettes.
- Multi channel pipette stand can hold 3 multi channel pipettes.

Pipette & Pipette Tips

